

Government of India
Directorate General of Health Services
Central Drugs Standard Control Organisation
Ministry of Health and Family welfare

FREE SALE CERTIFICATE

FSC Number : FSC/MD/2022/000123

Dated : 05/05/2022

This is to certify that M/s. Devon Innovations Private Limited situated at address No.27A, Near State Bank of India, Electronics city, Phase-1, Hosur road, Bangalore, Bengaluru (Bangalore) Urban Karnataka, 560100, 080-28522367, 368, India is holding Manufacturing License no. MFG/MD/2021/000159 issued on 09-JUN-2021 Under this license, the said firm is permitted to manufacture the products for their manufacture and sale in domestic market as per Indian Law and for export as per law of importing country (ies):-

Name of the products with their product name and dimensions if any mentioned below :

S.NO	Medical Device Details
1	1. Generic Name :Ureteral Indwelling Double Pigtail Stent/ Set (Short/Long Term) – with/without Hydrophilic coated(DJS)(OEO, BEO, OML OEO, OML BEO, BML BEO) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class C 5. Shelf life :5 Years 6. Sterile/Non-Sterile: Sterilized 7. Intended Use :Used for temporary internal drainage from the Ureteropelvic junction to the bladder 8. Material of Construction: Polyurethane/Carbothane 9. Dimension: 2.5F/6-32CM, 3F/6-32CM, 3.5F/6-32CM, 3.7/3.8F/6-32CM, 4F/8-32CM, 4.5F/10- 32CM, 4.7/4.8F/10-32CM, 5F/10-32CM, 5.5F/10-32CM, 6F/10-32CM, 6.5F/10-32CM, 7F/10-32CM, 7.5F/10-32CM, 8F/10-32CM, 9F/10-32CM, 10F/10-32CM and other dimensions 10. Model/Catalogue No./Name :NIL

VENUGOPA
L G SOMANI

Digitally signed by VENUGOPAL G SOMANI
DN: cn=IN, o=CENTRAL DRUGS STANDARD
CONTROL ORGANIZATION, ou=DRUGS
CONTROLLER GENERAL, (INDIA),
pseudoym=11204634621727ba3b2557487
32c456d49073e56953e9e9d5b8748ee3113
c, postalCode=110002, st=DELHI,
serialNumber=97241aa0c6ba41522026fca
0705074d4997e6a2f46dc8b81cf283adeaf5d,
cn=VENUGOPAL G SOMANI
Date: 2022.05.05 15:08:14 +05'30'

2	1. Generic Name :Ureteral Indwelling single Pigtail(Mono-J) Stent/ Set - with/without Hydrophilic coated 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class C 5. Shelf life :5 Years 6. Sterile/Non-Sterile: Sterilized 7. Intended Use :Used for temporary internal drainage from the Ureteropelvic junction to the bladder 8. Material of Construction: Polyurethane/Carbothane 9. Dimension: 2.5F/8-70CM, 3F/8-70CM, 3.5F/8-70CM, 3.7F/8-70CM, 4F/8-70CM, 4.5F/10-70CM, 4.7/4.8F/10-70CM, 5F/10-70CM, 5.5F/10-70CM, 6F/10-70CM, 6.5F/10-70CM, 7F/10-70CM, 8F/10- 70CM, 9F/10-70CM, 10F/10-70CM and other dimensions 10. Model/Catalogue No./Name :NIL
3	1. Generic Name :Biliary Stent/ Set- with/without Hydrophilic coated (BS)(BSAF,BSF, DPBS,BSHS,SPBS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class C 5. Shelf life :5 Years 6. Sterile/Non-Sterile: Sterilized 7. Intended Use :Used to drain obstructed biliary ducts 8. Material of Construction: Polyurethane 9. Dimension: 5F/3-15CM, 6F/3-15CM, 7F/3-15CM, 8F/3-15CM, 8.5F/3-15CM, 9F/3-15CM, 10F/3- 15CM and other dimensions 10. Model/Catalogue No./Name :NIL
4	1. Generic Name :Biliary Pancreatic Stent/Set-with/without Hydrophilic coated(BPS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class C 5. Shelf life :5 Years 6. Sterile/Non-Sterile: Sterilized 7. Intended Use :Used to drain obstructed biliary ducts 8. Material of Construction: Polyurethane 9. Dimension: 5F/3-15CM, 6F/3-15CM, 7F/3-15CM, 8F/3-15CM, 9F/3-15CM, 10F/3-15CM and other dimensions 10. Model/Catalogue No./Name :NIL
5	1. Generic Name :Ilyostomy stent with/without Hydrophilic coated(ILS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class C 5. Shelf life :5 Years 6. Sterile/Non-Sterile: Sterilized 7. Intended Use :Used for temporary drainage of urine from the ureteropelvic junction to the external urine collection bag during the condition of more infection in bladder and urethra 8. Material of Construction: Polyurethane/Carbothane 9. Dimension: 4F/60-70CM, 5F/60-70CM, 6F/60-70CM, 7F/60-70CM, 8F/60-70CM and other dimensions 10. Model/Catalogue No./Name :NIL



Government of India
Directorate General of Health Services
Central Drugs Standard Control Organisation
Ministry of Health and Family welfare

FREE SALE CERTIFICATE

FSC Number : FSC/MD/2022/000265

Dated : 08/09/2022

This is to certify that M/s. Devon Innovations Private Limited situated at address No.27A, Near State Bank of India, Electronics city, Phase-1, Hosur road, Bangalore ,Bengaluru (Bangalore) Urban Karnataka ,560100,080-28522367,368, India is holding Manufacturing License no. MFG/MD/2021/000297 issued on 13-SEP-2021 Under this license, the said firm is permitted to manufacture the products for their manufacture and sale in domestic market as per Indian Law and for export as per law of importing country (ies):-

Name of the products with their product name and dimensions if any mentioned below :

S.NO	Medical Device Details
1	1. Generic Name :Biliary Duct Drainage Catheter(BDDC) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for trans-hepatic biliary drainage. 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 6F/50CM,8F/50CM,10F/50CM, 12F/50CM,14F/50CM, 16F/50CM, 18F/50CM and other Dimensions 10. Model/Catalogue No./Name :NIL
2	1. Generic Name :Percutaneous Biliary Drainage Catheter/Set(PBDC) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Intended for percutaneous drainage in a variety of drainage applications 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 5F/32,35,40,42,45CM, 6F/32,35,40,42,45CM, 7F/32,35,40,42,45CM , 8F/32,35,40,42,45CM, 10F/32,35,40,42,45CM , 12F/32,35,40,42,45CM, 14F/32, 35,40,42,45CM, 16F/32,35,40,42,45CM, 18F/32,35,40,42,45CM and other Dimensions 10. Model/Catalogue No./Name :NIL

BHAGOJI TAYAPPA
KHANAPURE

Digitally signed by BHAGOJI
TAYAPPA KHANAPURE
Date: 2022.09.08 12:20:02
+05'30'

3	1. Generic Name :Tieman Catheter(TIEM) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to aid in urination for men who have enlarged prostates or suffer from benign prostatic hyperplasia 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 8F/20,40CM, 10F/20,30,40CM, 12F/20,30,40CM 14F/20,30,40CM, 16F/20,40CM, 18F/20,40CM, 20F/20,40CM, 22F/20,40CM, 24F/20,40CM and other Dimensions 10. Model/Catalogue No./Name :NIL
4	1. Generic Name :Endoureterotomy Stent/ Set(Internal)(EPIN) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :used for temporary drainage from the ureteropelvic junction to the bladder following incision of a stricture 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 8/4F-5/16-16CM, 10/5F-13/22,24,26CM, 10/6-13/26-22,24,26,28CM 10/5.5-13/26 CM 12/6-13/26-22,24,26,28CM 12/7-13/26-24,26CM 14/7-13/26-22,24,26,28 CM 14/6-13/26-22,24,26,28CM,16/8F-13/26-26cm and other Dimensions 10. Model/Catalogue No./Name :NIL
5	1. Generic Name :Internal External Drainage Catheter/ Set(IEDC) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for temporary drainage of the biliary duct 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 6F/40,45CM, 7F/30,40,45CM, 8F/40CM, 9F/40CM, 10F/40CM, 12F/40CM, 14F/40CM,16F/40CM, 18F/40CM, and other Dimensions 10. Model/Catalogue No./Name :NIL
6	1. Generic Name :Cutting Loop – (Single / Double Stem)-CLSS/CLDS 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :4 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to resect tissue and provide a smooth tissue cut 8. Material of Construction: Stainless steel 9. Dimension: 24F-30cm and other dimensions 10. Model/Catalogue No./Name :NIL

7	1. Generic Name :Multi – Band Ligator Barrel(LMB) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :4 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to endoscopically ligate esophageal varices at or above the gastroesophageal junction or to ligate internal hemorrhoids 8. Material of Construction: Latex free protein 9. Dimension: 6 bands-122 cm 10. Model/Catalogue No./Name :NIL
8	1. Generic Name :Foreign Body Retriever Dormia Basket(SB)-FLT/SBHEL/SEG/NIT/GRASPER 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :4 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for stone manipulation and retrieval / removal of Ureteral Stones from the Bladder, Renal Pelvis and Ureter 8. Material of Construction: Stainless steel 9. Dimension: 3F-70,90,110CM, 4F-70,90,110CM, 5F-70,90,110CM, 6F-70,90,110CM 10. Model/Catalogue No./Name :NIL
9	1. Generic Name :Guidewires/zebra Guidewires(GW) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3&5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Guidewires facilitate the delivery of wide variety of Catheters and Stents to a Medical procedure site within the body as a transient use 8. Material of Construction: Stainless steel & Nitinol 9. Dimension: 0.018IN/80,150CM, 0.021IN/80,150CM 0.025IN/80,150CM, 0.028IN/80,150CM 0.032IN/80,150,260,400CM, 0.035IN/80,150,260,400CM 0.038IN/80,150,260,400CM, and other Dimensions 10. Model/Catalogue No./Name :NIL
10	1. Generic Name :Vario Drain(VD) 2. Brand Name :VARIO DRAIN 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Vario Drain is a Closed Wound Suction Drainage system used after the surgery (Post Surgery) for removing the waste blood & puss from the closed wound. 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 7F-620mm,9F-620mm,12F-620mm,14F-620mm,16F-620mm,18F-620mm and other dimensions 10. Model/Catalogue No./Name :NIL

Competent Authority



Government of India
Directorate General of Health Services
Central Drugs Standard Control Organisation
Ministry of Health and Family welfare

FREE SALE CERTIFICATE

FSC Number : FSC/MD/2022/000266

Dated : 08/09/2022

This is to certify that M/s. Devon Innovations Private Limited situated at address No.27A, Near State Bank of India, Electronics city, Phase-1, Hosur road, Bangalore, Bengaluru (Bangalore) Urban Karnataka, 560100, 080-28522367, 368, India is holding Manufacturing License no. MFG/MD/2021/000297 issued on 13-SEP-2021 Under this license, the said firm is permitted to manufacture the products for their manufacture and sale in domestic market as per Indian Law and for export as per law of importing country (ies):-

Name of the products with their product name and dimensions if any mentioned below :

S.NO	Medical Device Details
1	1. Generic Name :Percutaneous Pigtail Nephrostomy catheter/Set(PTDC/ PCNS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for temporary or permanent drainage of urine from the kidney 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 5F/25,30,35CM, 6F/22,25,30,35CM, 7F/22,25,30,35CM, 8F/22, 25,30,35,40CM, 8.5F/25,30,35CM, 9F/22,25,30,35CM 10F/22,25,30,35CM, 12F/22, 25,30,35CM 14F/22,25,30,35CM, 16F/22,25,30,35CM 18F/25,30,35CM and other Dimensions 10. Model/Catalogue No./Name :NIL
2	1. Generic Name :IUI Catheter/Cannula(IUI) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for introduction of washed spermatozoa into the uterine cavity. Available in both straight and curved options. 8. Material of Construction: Medical Grade Polypropylene, low density polyethylene 9. Dimension: 6F-11CM, 6/8-11 cm/ 7.5 cm, 12 cm/ 9cm, 12 cm/ 15cm, 14 cm/ 10.5 cm, 17 cm/ 12 cm, 17.3 cm/ 13 cm and other dimensions 10. Model/Catalogue No./Name :NIL

3	1. Generic Name :Percutaneous Malecot Nephrostomy Catheter/Set(MNDC/PCMS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 & 5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for nephrostomy drainage in the renal pelvis. 8. Material of Construction: Polyvinyl chloride/ Ethylene-vinyl acetate & Medical Grade Polyurethane 9. Dimension: 8F/20,30CM, 10F/16,20,30CM, 12F/20,30CM 14F/20,30CM, 16F/20,30CM, 18F/20,30CM 20F/30,40CM, 22F/30CM, 24F/30CM 26F/30CM, 28F/30CM, 30F/30CM and other Dimensions 10. Model/Catalogue No./Name :NIL
4	1. Generic Name :Nephrostomy Catheter(NCD) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for drainage of the urine from the bladder at regular intervals 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 8F/30,40CM, 10F/30,35,40CM, 12F/20,30,40CM 14F/20,30,40CM, 16F/30,40CM, 18F/30,40CM 20F/30,40CM, 22F/30,40CM, 24F/30,40CM,26F/30,40CM, 28F/30,40CM 30F/30,40CM and other Dimensions 10. Model/Catalogue No./Name :NIL
5	1. Generic Name :Suprapubic Cystostomy Catheter/ Set(SPCS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to provide bladder drainage by percutaneous placement of Malecot catheter 8. Material of Construction: Medical Grade Ethylene-vinyl acetate 9. Dimension: 8F/20CM,10F/20CM, 12F/20,30CM, 14F/20CM, 16F/20CM,18F/20CM and other Dimensions 10. Model/Catalogue No./Name :NIL
6	1. Generic Name :Percutaneous Pigtail Suprapubic Catheter / Set(TDCS) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to provide bladder drainage by percutaneous placement of a pigtail catheter. 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 5F/12,22,24,26CM, 6F/12,22CM, 7F/12,22CM, 8F/12,22CM, 8.5F/12,22CM, 9F/12,22CM 10F/12,18,22CM, 12F/12,13,22CM, 14F/12,22CM, 15F/12,22CM, 16F/12,22CM, 18F/12,22CM, and other Dimensions 10. Model/Catalogue No./Name :NIL

7	1. Generic Name :Ureteral Catheter(UCB/ UOB (WT/AN/PT/CT/MH) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for drainage and retrograde pyelogram 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 3F/45,70CM, 4F/45,50,70CM, 5F/40,45,70CM 6F/45,70CM, 7F/45,70CM, 8F/45,70CM and other Dimensions 10. Model/Catalogue No./Name :NIL
8	1. Generic Name :Urinary Catheter(CIC-M/F) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for draining urine from the bladder in regular intervals. 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 6F/15,25,30,40CM, 8F/15,25,30,40CM, 10F/15,18,25,27,30,40,44CM, 12F/15,18,20,25,27,30,40,44CM, 14F/15,18,20,25,27,30,40,44CM 16F/15, 18,20,25,30,40,44CM, 18F/15,25,30,40CM, 20F/15,25,30,40CM 22F/15,25,30,40CM, 24F/15,25,30,40CM and other Dimensions 10. Model/Catalogue No./Name :NIL
9	1. Generic Name :Dilator/Set(UD/FDS/ FUD/FD/UDS/SD/MDF/MDM/MDP) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for dilation of the ureter prior to ureteroscopy and/or stone manipulation 8. Material of Construction: Medical Grade Polyurethane, Poly propylene, Ethylene-vinyl acetate, Poly propylene+Ethylene-vinyl acetate 9. Dimension: 5F/15,20,22,25,30,70CM, 6F/15,20,22,25,30,70CM, 7F/20,22,25,30,70CM, 8F/20,22,25,30,70CM, 9F/20,22,25,30,70CM, 10F/20,22,25,30,70CM, 12F/20, 22,25,30,70CM, 14F/20,22,25,30,70CM, 16F/20,22,25,30,70CM, 18F/20,22,25,30,70CM, 20F/20,22,25,30,70CM, 22F/20,22,25,30,70CM, 24F/20,22,25,30,70CM, 26F/20, 22,25,30,70CM 28F/20,22,25,30,70CM, 30F/20,22,25,30,70CM, (12-16)F/8cm,(14-24) F/9cm,(6-16)F/11cm and other Dimensions 10. Model/Catalogue No./Name :NIL

BHAGOJI TAYAPPA
KHANAPURE

Digitally signed by BHAGOJI
TAYAPPA KHANAPURE
Date: 2022.09.08 12:22:02 +05'30'

10	1. Generic Name :Amplatz Renal Introducer/ Set(ARI) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for nephrostomy tract dilation and the placement of a stricture 8. Material of Construction: Medical Grade Poly Tetra Fluoro Ethylene 9. Dimension: 8F/30CM,10F/30CM,12F/30CM,14F/30CM,16F/30CM,18F/30CM,20F/30CM,22F/30CM,24F/30CM,26F/30CM,28F/30CM,30F/30CM,32F/30CM,34F/30CM and other dimensions 10. Model/Catalogue No./Name :NIL
----	---



BHAGOJI TAYAPPA KHANAPURE Digitally signed by BHAGOJI TAYAPPA KHANAPURE
Date: 2022.09.08 12:22:59 +05'30'

Competent Authority

Government of India
Directorate General of Health Services
Central Drugs Standard Control Organisation
Ministry of Health and Family welfare

FREE SALE CERTIFICATE

FSC Number : FSC/MD/2023/000496

Dated : 26/09/2023

This is to certify that M/s. DEVON INNOVATIONS PVT LTD ,Devon Innovations Private Limited situated at address No.27A, Near State Bank of India, Electronics city, Phase-1,Hosur road, Bangalore ,Bengaluru (Bangalore) Urban Karnataka ,560100,080-28522367,368, India is holding Manufacturing License no. MFG/MD/2021/000297 issued on 13-SEP-2021 Under this license, the said firm is permitted to manufacture the products for their manufacture and sale in domestic market as per Indian Law and for export as per law of importing country (ies):-

Name of the products with their product name and dimensions if any mentioned below :

S.NO	Medical Device Details
1	1. Generic Name :External Biliary Drainage Catheter(EBDC) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for temporary drainage of the biliary duct 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 7F/25,50CM, 8F/25,40,50CM, 10F/25,40,50CM,12F/25,50CM,14F/25,50CM and other Dimensions 10. Model/Catalogue No./Name :NIL
2	1. Generic Name :Naso Biliary Drainage Catheter(NBDC) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for temporary drainage of the biliary duct 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 5F/200,220,250CM, 6F/200,220,250CM, 7F/200,220,250CM 8F/200,220,250CM, 10F/200,220,250CM and other Dimensions 10. Model/Catalogue No./Name :NIL

BHAGOJI TAYAPPA KHANAPURE

Digitally signed by BHAGOJI TAYAPPA KHANAPURE
Date: 2023.09.26 11:50:04 +05'30'

3	1. Generic Name :Urethral Stent/ Set(US)-(Short/Long Term) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :used for stenting the urethra after hypospadias or epispadias repair and to provide postoperative drainage of the bladder 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 5F/20,25CM, 6F/20,25CM, 7F/20,25CM, 8F/20,25CM, 10F/20,25CM, 12F/20,25CM, 14F/20,25CM, 16F/20,25CM and other dimensions 10. Model/Catalogue No./Name :NIL
4	1. Generic Name :Endopyelotomy Stent/ Set(EPEX)(External) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :used for stenting the Ureteropelvic junction and the ureter 8. Material of Construction: Medical Grade Polyurethane 9. Dimension: 8/4F-26/26--22,24,26CM, 10/5F-26/26-22,24,26CM, 10/6F-26/26-22,24,26CM, 12/6F-26/26-22,24,26CM, 14/7F-26/26-22,24,26CM, 16/7F-26/26-22,24,26CM, 5/14F-26/26-22,24,26CM, 6/14F-26/26-22,24,26CM and other Dimensions 10. Model/Catalogue No./Name :NIL
5	1. Generic Name :Endoscopic Jejunal Feeding Tube(EJFT) 2. Brand Name :DEVON 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 & 5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used through the channel of the endoscope to provide short-duration enteral access for delivery of nutrition or medication to the small bowel. 8. Material of Construction: Medical Grade Polyvinyl chloride & Medical Grade Polyurethane 9. Dimension: 10F/150CM, 12F/150CM, 14F/150CM, 16F/150CM and other Dimensions 10. Model/Catalogue No./Name :NIL
6	1. Generic Name :Klearsite Cannula(KLST) 2. Brand Name :KLEARSITE 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Klearsite is used during the surgery (Inter Operative) for removing the waste blood from the wound at the time of surgery & collect that into the suction apparatus. 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 10mm & 10.0mm/TC 10. Model/Catalogue No./Name :NIL

7	1. Generic Name :Biliary Duct Drainage Catheter(BDDC)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for trans-hepatic biliary drainage 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 5F/50CM,6F/50CM,7F/50CM,8F/50CM,10F/50CM, 12F/50CM,14F/50CM, 16F/50CM,18F/50CM and other Dimensions 10. Model/Catalogue No./Name :NIL
8	1. Generic Name :Connecting Tube(CT) 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for connecting Catheters to the urine bag 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 14F-10,30 cm, 16F-10 cm, 18F-10cm, 24F-10 cm and other Dimensions 10. Model/Catalogue No./Name :NIL
9	1. Generic Name :Feeding Tube(FT) 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to supply nutrition 8. Material of Construction: Medical Grade Ethylene-vinyl acetate 9. Dimension: 10F/125 cm,12F/125 cm, 14F/125 cm,16F/125 cm and other Dimensions 10. Model/Catalogue No./Name :NIL
10	1. Generic Name :Re-entry Malecot catheter(REMC) 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for internal as well as external drainage with combination of nephrostomy catheter with stent 8. Material of Construction: Medical Grade Ethylene-vinyl acetate 9. Dimension: 10/5F-22/26 cm,24/26 cm ,26 /26 cm, 10/6F-22/26cm,24/26cm,26 /26cm, 12/6F- 22/26cm,24/26cm,26 /26cm,30/17cm, 14/6F-22/26cm,24/26cm, 26 /26cm,30/17cm, 14/7F- 22/26cm,24/26cm, 26 /26cm,30/17cm, 16/7F-22/26cm, 24/26cm, 26 /26cm,30/17cm, 16/8F- 22/26cm,24/26cm, 26 /26cm,30/17cm, 18/8F- 22/26cm,24/26cm, 26 /26cm,30/17cm, 20/8F- 22/26cm,24/26cm, 26 /26cm,30/17cm and other Dimensions 10. Model/Catalogue No./Name :NIL

Government of India
Directorate General of Health Services
Central Drugs Standard Control Organisation
Ministry of Health and Family welfare

FREE SALE CERTIFICATE

FSC Number : FSC/MD/2022/000144

Dated : 25/05/2022

This is to certify that M/s. Devon Innovations Private Limited situated at address No.27A, Near State Bank of India, Electronics city, Phase-1, Hosur road, Bangalore ,Bengaluru (Bangalore) Urban Karnataka ,560100,080-28522367,368, India is holding Manufacturing License no. MFG/MD/2021/000297 issued on 13-SEP-2021 Under this license, the said firm is permitted to manufacture the products for their manufacture and sale in domestic market as per Indian Law and for export as per law of importing country (ies):-

Name of the products with their product name and dimensions if any mentioned below :

S.NO	Medical Device Details
1	1. Generic Name :Percutaneous Pigtail Nephrostomy catheter/Set(PTDC/ PCNS)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for temporary or permanent drainage of urine from the kidney 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 5F/25,30,35,40CM, 6F/25,30,35,40CM, 7F/25,30,35,40CM, 8F/25,30,35,40CM, 8.5F/25,30,35,40CM, 9F/25,30,35,40CM 10F/25,30,35,40CM, 12F/25,30,35,40CM 14F/25,30,35,40CM, 16F/25,30,35,40CM 18F/25,30,35,40 CM and other Dimensions 10. Model/Catalogue No./Name :NIL

BHAGOJI TAYAPPA
KHANAPURE

Digitally signed by BHAGOJI
TAYAPPA KHANAPURE
Date: 2022.05.25 18:21:25 +05'30'

2	1. Generic Name :Percutaneous Malecot Nephrostomy Catheter/Set(MNDC)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 & 5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for nephrostomy drainage in the renal pelvis 8. Material of Construction: Polyvinyl chloride/ Ethylene-vinyl acetate & Medical Grade Polyurethane 9. Dimension: 8F/20,25,30CM, 10F/20,25,30CM, 12F/20,25,30CM 14F/20,25,30CM, 16F/20,25,30CM, 18F/20,25,30CM 20F/25,30CM, 22F/25,30CM, 24F/25,30CM 26F/25, 30CM, 28F/25,30CM, 30F/25,30CM and other Dimensions 10. Model/Catalogue No./Name :NIL
3	1. Generic Name :Guidewires(GW)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Guidewires facilitate the delivery of wide variety of Catheters and Stents to a Medical procedure site within the body as a transient use 8. Material of Construction: Stainless Steel, Nitinol-Hydrophilic coating 9. Dimension: 0.018IN/80,150CM, 0.021IN/80,150CM 0.025IN/80,150CM, 0.028IN/80, 150CM 0.032 IN/80,150,260,400CM, 0.035IN/80,150,260,400CM 0.038IN/80, 150,260,400CM, and other Dimensions 10. Model/Catalogue No./Name :NIL
4	1. Generic Name :Ureteral Access sheath(UAS)- Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to establish a conduit during endoscopic urological procedures facilitating the passage of endoscopes and other instruments into the urinary tract 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating,Stainless Steel 9. Dimension: 9.5Fr/11.5Fr – 35,38,45,55 cms, 10.7Fr/12.5Fr – 35,38,45,55 cms, 12Fr/14Fr – 35,38,45,55 cms and other Dimensions 10. Model/Catalogue No./Name :NIL
5	1. Generic Name :TUR Set(TUR) 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :TUR set used for Endoscopic irrigation during Trans Urethral Resection of prostate gland 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 22F-143 CM, 24F-143 CM and other Dimensions 10. Model/Catalogue No./Name :NIL

6	1. Generic Name :Abdominal Drainage Kit(ADK) 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Specially designed for post operative abdominal drainage 8. Material of Construction: Medical Grade Polyvinyl chloride 9. Dimension: 18F-22 cm, 20F-22, 22F-22 cm, 24F-22 cm,28F-22 cm, 32F-22 cm and other Dimensions 10. Model/Catalogue No./Name :NIL
7	1. Generic Name :Percutaneous Nephrolithotomy Set(PCNLS) 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :These are convenience kits containing multiple devices to perform removal of stones in the kidney 8. Material of Construction: Medical Grade Ethylene-vinyl acetate+ Polypropylene 9. Dimension: 8F-20,22 cm 10F-20,22 cm 12F-20,22 cm 14F-20,22 cm 16F-20,22 cm 18F-20,22 cm and other Dimensions 10. Model/Catalogue No./Name :NIL



BHAGOJI TAYAPPA
KHANAPURE

Competent Authority

Digitally signed by BHAGOJI
TAYAPPA KHANAPURE
Date: 2022.05.25 18:22:09 +05'30'

Government of India
 Directorate General of Health Services
 Central Drugs Standard Control Organisation
 Ministry of Health and Family welfare

FREE SALE CERTIFICATE

FSC Number : FSC/MD/2023/000495

Dated : 26/09/2023

This is to certify that M/s. DEVON INNOVATIONS PVT LTD ,Devon Innovations Private Limited situated at address No.27A, Near State Bank of India, Electronics city, Phase-1,Hosur road, Bangalore ,Bengaluru (Bangalore) Urban Karnataka ,560100,080-28522367,368, India is holding Manufacturing License no. MFG/MD/2021/000297 issued on 13-SEP-2021 Under this license, the said firm is permitted to manufacture the products for their manufacture and sale in domestic market as per Indian Law and for export as per law of importing country (ies):-

Name of the products with their product name and dimensions if any mentioned below :

S.NO	Medical Device Details
1	1. Generic Name :Percutaneous Biliary Drainage Catheter(PBDC)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Intended for percutaneous drainage in a variety of drainage applications 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 5F/32,35,40,42,45CM, 6F/32,35,40,42,45CM, 7F/32,35,40,42,45CM , 8F/32, 35,40,42,45CM, 10F/32,35,40,42,45CM , 12F/32,35,40,42,45CM, 14F/32, 35,40,42,45CM, 16F/32, 35,40,42,45CM 18F/32,35,40,42,45CM and other Dimensions 10. Model/Catalogue No./Name :NIL
2	1. Generic Name :Suprapubic Cystostomy Catheter/ Set(SPCS)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to provide bladder drainage by percutaneous placement of Malecot catheter 8. Material of Construction: Medical Grade Ethylene-vinyl acetate-Hydrophilic coating 9. Dimension: 8F/20CM,10F/20CM, 12F/20CM 14F/20CM, 16F/20CM,18F/20CM and other Dimensions 10. Model/Catalogue No./Name :NIL

3	1. Generic Name :Percutaneous Pigtail Suprapubic Catheter / Set(TDCS)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used to provide bladder drainage by percutaneous placement of a pigtail catheter 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 5F/12,22,24,26CM,6F/12,22CM, 7F/12,22CM, 8F/12,22,25CM, 8.5F/12, 22CM, 9F/12,22CM 10F/12,22CM, 12F/12,22CM, 14F/12,22CM, 16F/12,22CM, 18F/12,22CM and other Dimensions 10. Model/Catalogue No./Name :NIL
4	1. Generic Name :Ureteral Catheter(UCB/ UOB-WT/AN/PT/CT/MH)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for drainage and retrograde pyelogram 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 3F/45,70CM, 4F/45,70CM, 5F/45,70CM 6F/45,70CM, 7F/45,70CM, 8F/45, 70CM and other Dimensions 10. Model/Catalogue No./Name :NIL
5	1. Generic Name :Urinary Catheter(CIC-M/F)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for draining urine from the bladder in regular intervals 8. Material of Construction: Medical Grade Polyvinyl chloride-Hydrophilic coating 9. Dimension: 6F/15,25,30,40CM,8F/15,25,30,40CM, 10F/15,25,30,40CM, 12F/15, 25,30,40CM, 14F/15,25,30,40CM 16F/15,25,30,40CM, 18F/15,25,30,40CM, 20F/15, 25,30,40CM 22F/15,25,30,40CM, 24F/15,25,30,40 CM and other Dimensions 10. Model/Catalogue No./Name :NIL

BHAGOJI TAYAPPA KHANAPURE Digitally signed by BHAGOJI TAYAPPA KHANAPURE
Date: 2023.09.26 11:49:08 +05'30'

6	1. Generic Name :Dilator/Set(UD/FDS/ FUD/FD/UDS/SD)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :3 years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for dilation of the ureter prior to ureteroscopy and/or stone manipulation 8. Material of Construction: Medical Grade Poly propylene, Ethylene-vinyl acetate,Poly propylene+Ethylene-vinyl acetate,Polyurethane-Hydrophilic coating 9. Dimension: 5F/20,22,25,30,70CM, 6F/20,22,25,30,70CM, 7F/20,22,25,30,70CM, 8F/20, 22,25,30,70CM, 9F/20,22,25,30,70CM, 10F/20,22,25,30,70CM, 12F/20, 22,25,30,70CM, 14F/20, 22,25,30,70CM, 16F/20,22,25,30,70CM, 18F/20,22,25,30,70CM, 20F/20,22,25,30,70CM, 22F/20, 22,25,30,70CM, 24F/20,22,25,30,70CM, 26F/20, 22,25,30,70CM 28F/20,22,25,30,70CM, 30F/20, 22,25,30,70CM and other Dimensions 10. Model/Catalogue No./Name :NIL
7	1. Generic Name :Urethral Stent/ Set(US)- (Short/Long Term) Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :used for stenting the urethra after Hypospadias or Epispadias repair and to provide post operative drainage of the bladder 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 5F/20,25CM, 6F/20,25CM, 7F/20,25CM, 8F/20,25CM, 10F/20,25CM, 12F/20,25CM, 14F/20,25CM, 16F/20,25CM and other dimensions 10. Model/Catalogue No./Name :NIL
8	1. Generic Name :Endopyelotomy Stent/ Set(EPEX)(External)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :used for stenting the Ureteropelvic junction and the ureter 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 10/5F-26/26,24/26CM, 12/6F-26/26,24/26CM, 14/7F-26/26,24/26CM, 5/14F- 26/26,24/26CM, 6/14F-26/26,24/26CM and other Dimensions 10. Model/Catalogue No./Name :NIL
9	1. Generic Name :Endoureterotomy Stent/ Set(Internal)(EPIN)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :used for temporary drainage from the ureteropelvic junction to the bladder following incision of a stricture 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 8/4F-5/16, 10/5F-13/22,24,26CM, 10/6-13/22,24,26,28CM 10/5.5-13/26 CM 12/6-13/22,24,26,28CM 12/7-13/24,26CM 14/7-13/22,24,26,28CM 14/6-13/ 22,24,26,28CM and other Dimensions 10. Model/Catalogue No./Name :NIL

10	1. Generic Name :Internal External Drainage Catheter/ Set(IEDC)-Hydrophilic coated 2. Brand Name :Devon 3. Notified Category :Catheters 4. Class of Medical Device :Class B 5. Shelf life :5 Years 6. Sterile/Non-Sterile:Both 7. Intended Use :Used for temporary drainage of the biliary duct 8. Material of Construction: Medical Grade Polyurethane-Hydrophilic coating 9. Dimension: 6F/40CM, 7F/40CM, 8F/40CM, 9F/40CM, 10F/40CM, 12F/40CM, 14F/40CM, 16F/40CM, 18F/40CM, 20F/40CMand other Dimensions 10. Model/Catalogue No./Name :NIL
----	---

BHAGOJI TAYAPPA KHANAPURE

Digitally signed by BHAGOJI TAYAPPA KHANAPURE
Date: 2023.09.26 11:49:38 +05'30'

Competent Authority